



AIROS 6 Sequential Compression Device

Reliable, Easy-to-Use 6-Chamber Device Ideal for the Medicare Population

The FDA-cleared AIROS® 6 Sequential Compression Device delivers compression therapy to help patients manage lymphatic and venous disorders. The device provides effective intermittent pneumatic compression (IPC) therapy for treatment of the lower and upper extremities.

The AIROS 6 provides gradient sequential compression that increases blood flow and encourages extracellular fluid clearance. The device provides cycles of compressed air at certain adjustable pressures and then inflates and deflates garments to move fluid.

Overview & Features

Regulatory Clearances: FDA 510(k) & Medicare (HCPCS Code E0651); also reimbursable through commercial insurance

Usage Tracking: Built-in timer records total usage time to measure patient compliance

Easy Operation: Large, clear LCD screen and one-button operation for starting and stopping therapy

Physical Benefits: Reduce pain and swelling, increase blood flow, move lymphatic fluid

Various Sizes of Comfortable Garments: 10 garments for lower extremity treatment & 6 garments for upper extremity treatment, including arm/shoulder/chest (NOTE: Not available for use with double extra wide leg or pants garments. These garments must be used with the AIROS 6P device)

Lymphedema Preparation Treatment (LPT) Modes: Optional Clearing Mode (mimics lymphatic drainage) or Focal Mode (fibrotic tissue relief) can be used before main Gradient Mode therapy

Powerful Vacuum: Air is quickly vacuumed from the garment after therapy is completed to assist in garment removal



AIROS 6/6P
Instructions for
Use Guide (PDF)



Additional Features for AIROS 6P Model

✓ **Universal Electric:** The AIROS 6P Model features 100-240 VAC, 50/60Hz electrical rating, which allows the device to be used anywhere in the world.

✓ **Peristaltic Mode:** The Peristaltic Mode provides a wave-like feeling to the patient, as the pressure is not held in all chambers throughout a cycle of treatment. Some lymphedema therapists prefer this type of therapy for patients with upper extremity swelling.

The user can set the starting pressure in Chamber 1 near the hand. Pressure decreases at a fixed 5% gradient distally to proximally. Chambers 1 through 3 first inflate to their target pressures sequentially and hold. Then Chamber 1 deflates, relieving pressure near the hand and wrist, and Chamber 4 inflates to its target pressure and holds. This process continues as the 6 chambers are inflated moving up the extremity.

NOTE: The AIROS 6P must be used for the pants garments which are prescribed for truncal swelling, and some of the larger sized leg garments.

AIROS 6P Sequential Compression Device

6-Chamber E0651 Model for Use with Peristaltic Mode, Truncal Treatment & Larger Leg Garments

The AIROS 6P model is based on the AIROS 6 device and provides all the same therapeutic benefits and most of the same technical characteristics. The main difference is that the 6P features an additional primary operating mode, the Peristaltic Mode, for use with the Arm and Arm plus garments (upper extremity only).



AIROS 6/6P Garment Use Reference Guide

Garment Type/Model	AIROS 6	AIROS 6P
LOWER EXTREMITY		
6-SL01, 6-SWL01, 6-SXWL01	YES	YES
6-SDXWL01	NO	YES
6-ML01, 6-MWL01, 6-MXWL01	YES	YES
6-MDXWL01	NO	YES
6-LL01, 6-LWL01	YES	YES
6-LXWL01	NO	YES
6-HL01, 6-HWL01	YES	YES
UPPER EXTREMITY		
6-MA01, 6-LA01	YES, GRADIENT MODE ONLY	YES, GRADIENT OR PERISTALTIC MODE
6-APL01, 6-APLP01, 6-APR01, 6-APRP01	YES, GRADIENT MODE ONLY	YES, GRADIENT OR PERISTALTIC MODE
TRUNCAL (PANTS)		
6-SP01, 6-MP01, 6-LP01	NO	YES

AIROS 6/6P Specifications & Information Comparison

Model	AIROS 6	AIROS 6P
Dimensions	10.2"W x 10.5" D x 4.3"H	10.2"W x 10.5" D x 4.3"H
Weight	6.4 lbs.	6.4 lbs.
Electrical	120 VAC, 60Hz, 100VA MAX (USA Only)	100-240 VAC, 50/60Hz, 100VA MAX (USA & International)
Primary Operating Modes	1 (Gradient Mode)	2 (Gradient Mode and Peristaltic Mode)
Secondary Lymphedema Preparation Treatment (LPT) Modes	2 (Focal Mode and Clearing Mode)	2 (Focal Mode and Clearing Mode)
Pressure Range Setting (Primary Operating Modes)	30-80mmHg (increments of 1mmHg)	30-80mmHg (increments of 1mmHg)
Cycle Time (Inflation) Setting	Gradient Mode and Pressure Mode: 40 or 60 seconds	Gradient Mode and Pressure Mode: 40 or 60 seconds Peristaltic Mode: No setting
Therapy Time Setting	15, 30, 45 or 60 minutes	15, 30, 45 or 60 minutes
Device Warranty	3 years	3 years
Error Indicators	3	3
Air Blocker (Unilateral Therapy)	1 included	1 included