

AIROS™ 6 SEQUENTIAL COMPRESSION DEVICE



Reliable, Intuitive 6-Chamber Device
Serving the Medicare Population

The FDA-cleared AIROS 6 Sequential Compression Device is a compression therapy device utilized for the treatment and management of lymphedema and venous insufficiencies. The device provides effective intermittent pneumatic compression (IPC) therapy for treatment of the legs and arms.

The AIROS 6 provides gradient sequential compression that increases blood flow and encourages extracellular fluid clearance. The device provides cycles of compressed air at certain adjustable pressures and then inflates and deflates garments to move fluid.

KEY FEATURES

- **Billing:** Medicare patients (Code E0651) or private insurance
- **Ease of Use Promoting Patient Compliance:** Large, clear LCD screen and one-button operation for starting and stopping therapy
- **Physical Benefits:** Reduce pain and swelling, increase blood flow, move lymphatic fluid
- **Various Sizes of Comfortable Garments:** 8 full leg sizes, 2 half leg sizes, and 2 arm sizes available to comfortably fit patients (all utilize 6 chambers)
- **Lymphedema Preparation Treatment (LPT) Modes:** Optional Clearing Mode (mimics lymphatic drainage) or Focal Mode (fibrotic tissue relief) can be used before main Gradient Mode therapy



SPECIFICATIONS

| | |
|----------------------------------|---------------------------------|
| Dimensions | 10.2" W x 10.5" D x 4.3" H |
| Weight | 6.4 lbs. |
| Primary Operating Mode | Gradient |
| Secondary Operating Modes (LPT) | Clearing or Focal |
| Pressure Range (Gradient Mode) | 30-80mmHg (increments of 1mmHg) |
| Cycle Time (Inflation) | 40 or 60 seconds |
| Deflation Time | 12 seconds |
| Therapy Time | 15, 30, 45 or 60 minutes |
| Warranty (Device) | 3 years |
| Alarms (Error Messages) | 3 |
| Air Blocker (Unilateral Therapy) | Included (1) |
| Electrical | 120 VAC, 60Hz, 100VA MAX |