

AIROS® 8-CHAMBER COMPRESSION GARMENTS



Comfortable, Reliable Garments for
Upper & Lower Extremity Therapy

The AIROS® 8-Chamber Compression Garments are intended for use with the FDA-cleared AIROS 8 Sequential Compression Device. Together, this medical technology provides compression therapy for patients living with lymphedema and venous insufficiencies.

The compression garments were designed with patient comfort and clinical effectiveness in mind. From durable material to precise sizing, our garments were carefully engineered to support patients with various needs.

KEY FEATURES

- **Various Sizes:** 8 full leg sizes and 2 arm sizes available to comfortably fit patients (all utilize 8 individual chambers)
- **Ease of Use:** Pull handles and O-ring zippers to facilitate easy application and removal
- **High-Quality Material:** Durable, easy-to-clean material also allows for simple cleaning
- **Reimbursement:** Commercial insurance & Medicare (codes below)
 - E0667 - Leg Garment
 - E0668 - Arm Garment



AIROS® 8-CHAMBER COMPRESSION GARMENTS



Comfortable, Reliable Garments for
Upper & Lower Extremity Therapy

AIROS® 8-Chamber Compression Garments Sizing Chart

LEG GARMENTS

Garment	Model	Inseam	Circumference (Ankle)	Circumference (knee)	Circumference (Thigh)
Small Leg	8-SL01	66cm (26")	53cm (21")	64cm (25")	76cm (30")
Small Wide Leg	8-SWL01	66cm (26")	61cm (24")	74cm (29")	89cm (35")
Small Extra Wide Leg	8-SXWL01	66cm (26")	71cm (28")	89cm (35")	102cm (40")
Medium Leg	8-ML01	81cm (32")	56cm (22")	69cm (27")	79cm (31")
Medium Wide Leg	8-MWL01	81cm (32")	64cm (25")	79cm (31")	89cm (35")
Medium Extra Wide Leg	8-MXWL01	81cm (32")	74cm (29")	91cm (36")	102cm (40")
Large Leg	8-LL01	91cm (36")	53cm (21")	69cm (27")	79cm (31")
Large Wide Leg	8-LWL01	91cm (36")	56cm (22")	79cm (31")	102cm (40")

ARM GARMENTS

Garment	Model	Inseam	Circumference (Axilla)	Circumference (Elbow)	Circumference (Wrist)
Medium Arm	8-MA01	74cm (29")	56cm (22")	48cm (19")	33cm (13")
Large Arm	8-LA01	81cm (32")	58cm (23")	51cm (20")	36cm (14")