

# Lower Extremity Garments

## Durable, Comfortable Garments for Effective Therapy

AIROS Medical lower extremity compression garments are available in both 6-chamber and 8-chamber configurations for use with AIROS 6/6P and AIROS 8/8P devices.

These compression garments were designed with patient comfort, ease of use, and clinical effectiveness in mind. From durable material to a wide variety of sizes, our garments were carefully engineered to support patients with various needs.

Models span small, medium, and large inseam lengths and wide sizes are also available. There are also 2 half leg garments for use only with the AIROS 6/6P device.



### Key Features

#### Various Sizes:

##### 6-Chamber Garments:

11 Full Leg sizes, 2 Half Leg sizes

##### 8-Chamber Garments:

11 Full Leg sizes

**Ease of Use:** Pull handles and O-ring zippers to facilitate easy application and removal

#### High-Quality Material:

Durable material and allows for simple cleaning

#### Reimbursement:

Commercial Insurance

Medicare billing codes are E0667 (Full Leg Garment) and E0669 (Half Leg Garment)



Lower Extremity  
Compression  
Garments

*NOTE: Please refer to the sizing chart as some garment models are only for use with the AIROS 6P and AIROS 8P devices.*





## AIROS 6-Chamber & 8-Chamber Lower Extremity Garment Sizing Chart

Garment	Model	Inseam	Circumference (Ankle)	Circumference (Knee)	Circumference (Thigh)
Small Leg	6-SL01/8-SL01	66cm (26")	53cm (21")	64cm (25")	76cm (30")
Small Wide Leg	6-SWL01/8-SWL01	66cm (26")	61cm (24")	74cm (29")	89cm (35")
Small Extra Wide Leg	6-SXWL01/8-SXWL01	66cm (26")	71cm (28")	89cm (35")	102cm (40")
Small Double Extra Wide Leg	6-SDXWL01/8-SDXWL01	66cm (26")	81cm (32")	107cm (42")	124cm (49")
Medium Leg	6-ML01/8-ML01	81cm (32")	56cm (22")	69cm (27")	79cm (31")
Medium Wide Leg	6-MWL01/8-MWL01	81cm (32")	64cm (25")	79cm (31")	89cm (35")
Medium Extra Wide Leg	6-MXWL01/8-MXWL01	81cm (32")	74cm (29")	91cm (36")	102cm (40")
Medium Double Extra Wide Leg	6-MDXWL01/8-MDXWL01	81cm (32")	81cm (32")	109cm (43")	124cm (49")
Large Leg	6-LL01/8-LL01	91cm (36")	53cm (21")	69cm (27")	79cm (31")
Large Wide Leg	6-LWL01/8-LWL01	91cm (36")	56cm (22")	79cm (31")	102cm (40")
Large Extra Wide Leg	6-LXWL01/8-LXWL01	91cm (36")	84cm (33")	109cm (43")	127cm (50")
Half Leg	6-HL01	46cm (18")	56cm (22")	71cm (28")	-
Half Wide Leg	6-HWL01	46cm (18")	61cm (24")	76cm (30")	-