

Lower Extremity Garments

Durable, Comfortable Garments for Effective Therapy

AIROS Medical lower extremity compression garments are available in both 6-chamber and 8-chamber configurations for use with AIROS 6/6P and AIROS 8/8P devices.

These compression garments were designed with patient comfort, ease of use, and clinical effectiveness in mind. From durable material to a wide variety of sizes, our garments were carefully engineered to support patients with various needs.

Models span small, medium, and large inseam lengths and wide sizes are also available. There are also 2 half leg garments for use only with the AIROS 6/6P device.



Key Features

Various Sizes:

6-Chamber Garments:
11 Full Leg sizes, 2 Half Leg sizes

8-Chamber Garments:
11 Full Leg sizes

Ease of Use: Pull handles and O-ring zippers to facilitate easy application and removal

High-Quality Material:
Durable material and allows for simple cleaning

Reimbursement:

Commercial Insurance
Medicare billing codes are E0667
(Full Leg Garment) and E0669
(Half Leg Garment)



Lower Extremity
Compression
Garments

NOTE: Please refer to the sizing chart as some garment models are only for use with the AIROS 6P and AIROS 8P devices.





AIROS
6-Chamber
Lower
Extremity
Garment

AIROS
8-Chamber
Lower
Extremity
Garment

**AIROS 6-Chamber & 8-Chamber
Lower Extremity Garment Sizing Chart**

Garment	Model	Inseam	Circumference (Ankle)	Circumference (Knee)	Circumference (Thigh)
Small Leg	6-SL/8-SL	66cm (26")	53cm (21")	64cm (25")	76cm (30")
Small Wide Leg	6-SWL/8-SWL	66cm (26")	61cm (24")	74cm (29")	89cm (35")
Small Extra Wide Leg	6-SXWL/8-SXWL	66cm (26")	71cm (28")	89cm (35")	102cm (40")
Small Double Extra Wide Leg	6-SDXWL/8-SDXWL	66cm (26")	81cm (32")	107cm (42")	124cm (49")
Medium Leg	6-ML/8-ML	81cm (32")	56cm (22")	69cm (27")	79cm (31")
Medium Wide Leg	6-MWL/8-MWL	81cm (32")	64cm (25")	79cm (31")	89cm (35")
Medium Extra Wide Leg	6-MXWL/8-MXWL	81cm (32")	74cm (29")	91cm (36")	102cm (40")
Medium Double Extra Wide Leg	6-MDXWL/8-MDXWL	81cm (32")	81cm (32")	109cm (43")	124cm (49")
Large Leg	6-LL/8-LL	91cm (36")	53cm (21")	69cm (27")	79cm (31")
Large Wide Leg	6-LWL/8-LWL	91cm (36")	56cm (22")	79cm (31")	102cm (40")
Large Extra Wide Leg	6-LXWL/8-LXWL	91cm (36")	84cm (33")	109cm (43")	127cm (50")
Half Leg	6-HL	46cm (18")	56cm (22")	71cm (28")	-
Half Wide Leg	6-HWL	46cm (18")	61cm (24")	76cm (30")	-