


Truncal Garments

Pants Configuration Features Bilateral Leg Treatment with Compression Chambers Above the Waist to Treat Abdominal Swelling

AIROS Medical's pants garments treat the entire lower body, delivering compression therapy to the legs, abdomen, lower back, and hip area. Lymphatic drainage in these areas is critical for many patients.

These garments are made from the same durable, comfortable material as our leg, arm and Arm Plus garments, and feature pull handles for easy application and removal. The pants garments are available in both 6-chamber (used with the AIROS 6P) and 8-chamber (used with the AIROS 8P) configurations.



Key Features

Full Lower Body Therapy:

Compression chambers deliver treatment to the legs, abdomen, lower back, and hip area.

Ease of Use: Pull handles and O-ring zippers allow easy application and removal.

High-Quality Material: Durable, comfortable material that is easy to clean.

Reimbursement: Medicare and commercial insurance options.



Truncal
Garments



AIROS 6-Chamber & 8-Chamber Truncal Garment Sizing Chart

Garment	Model	Device Compatibility	Inseam Length	Abdomen Length	Total Length	Circumference (Ankle)	Circumference (Knee)	Circumference (Thigh)	Circumference (Abdomen)
Small Pants (6-Chamber)	6-SP01	AIROS 6P	66 cm (26")	46 cm (18")	117cm (46")	53cm (21")	79 cm (31")	79 cm (31")	123 cm (48")
Medium Pants (6-Chamber)	6-MP01	AIROS 6P	81cm (32")	46 cm (18")	127cm (50")	64 cm (25")	81cm (32")	81cm (32")	125 cm (49")
Large Pants (6-Chamber)	6-LP01	AIROS 6P	91 cm (36")	46 cm (18")	142cm (56")	71 cm (28")	94 cm (37")	94 cm (37")	145 cm (57")
Small Pants (8-Chamber)	8-SP01	AIROS 8P	66 cm (26")	46 cm (18")	117cm (46")	53cm (21")	79 cm (31")	79 cm (31")	123 cm (48")
Medium Pants (8-Chamber)	8-MP01	AIROS 8P	81cm (32")	46 cm (18")	127cm (50")	64 cm (25")	81cm (32")	81cm (32")	125 cm (49")
Large Pants (8-Chamber)	8-LP01	AIROS 8P	91 cm (36")	46 cm (18")	142cm (56")	71 cm (28")	94 cm (37")	94 cm (37")	145 cm (57")