


Upper Extremity Garments

Advanced, Easy-to-Use Garments for Patients with Lymphedema in the Arm, Shoulder, and Chest

AIROS Medical's upper extremity garments provide an effective treatment option for patients with lymphedema in the upper extremities, especially for post-mastectomy/breast cancer patients with lymphedema.

These compression garments were designed so that patients themselves - even elderly ones - can put on the garments and remove them without the help of a family member or caregiver. This ease-of-use increases patient compliance and drives positive outcomes.



Key Features

Various Sizes

Arm Garments: Medium and Large inseam lengths

Arm Plus Garments: Available in left and right configurations, and in two colors

Both garment types available in 6-chamber (E0651) and 8-chamber (E0652) sizes

6-chamber garments can be used with the AIROS 6 (Gradient Mode) and the AIROS 6P (Gradient or Peristaltic Mode)

8-chamber garments can be used with the AIROS 8 (Gradient Mode or Pressure Mode) and the AIROS 8P (Gradient, Pressure, or Peristaltic Mode)

High-Quality Material: Durable material and allows for simple cleaning

Ease of Use:

Arm Plus garment features support padding in the back to aid in donning and removal

Arm garment simply slides on and reaches the axilla area

Reimbursement:

Commercial Insurance

Medicare billing code used for both Arm and Arm Plus garment is E0668

Quality Therapy:

Patients prescribed an Arm Plus garment benefit from quality therapy immediately, so no lengthy trials are needed for billing purposes



AIROS 6-Chamber & 8-Chamber
Upper Extremity Compression
Garments



AIROS Arm Plus Compression Garment

AIROS Arm Compression Garment

AIROS 6-Chamber & 8-Chamber Upper Extremity Garment Sizing Chart

Garment	Model	Inseam	Circumference (Body)	Circumference (Axilla)	Circumference (Elbow)	Circumference (Wrist)
Arm Plus Garments						
Arm Plus Left	6-APL/8-APL	74cm (29")	76-173cm (30"-68")	56cm (22")	51cm (20")	36cm (14")
Arm Plus Pink Left	6-APLP/8-APLP	74cm (29")	76-173cm (30"-68")	56cm (22")	51cm (20")	36cm (14")
Arm Plus Right	6-APR/8-APR	74cm (29")	76-173cm (30"-68")	56cm (22")	51cm (20")	36cm (14")
Arm Plus Pink Right	6-APRP/8-APRP	74cm (29")	76-173cm (30"-68")	56cm (22")	51cm (20")	36cm (14")
Arm Garments						
Medium Arm	6-MA/8-MA	74cm (29")	--	56cm (22")	48cm (19")	33cm (13")
Large Arm	6-LA/8-LA	81 cm (32")	--	58cm (23")	51cm (20")	36cm (14")

Our Arm and Arm Plus models can be used with the AIROS 6/6P or AIROS 8/8P devices. Patients utilizing the AIROS 6P or AIROS 8P device can utilize the Arm and Arm Plus garments with the Peristaltic Mode, if the patient's doctor or lymphedema therapist believe this wave-like pressure inflation would be better than traditional gradient sequential compression therapy.