


Truncal Garments

Pants Configuration Features Bilateral Leg Treatment with Compression Chambers Above the Waist to Treat Abdominal Swelling

AIROS Medical's pants garments treat the entire lower body, delivering compression therapy to the legs, abdomen, lower back, and hip area. Lymphatic drainage in these areas is critical for many patients.

These garments are made from the same durable, comfortable material as our leg, arm and Arm Plus garments, and feature pull handles for easy application and removal. The pants garments are available in both 6-chamber (used with the AIROS 6P) and 8-chamber (used with the AIROS 8P) configurations.



Key Features

Full Lower Body Therapy:

Compression chambers deliver treatment to the legs, abdomen, lower back, and hip area.

Ease of Use:

Pull handles and O-ring zippers allow easy application and removal.

High-Quality Material:

Durable, comfortable material that is easy to clean.

Reimbursement:

Medicare and commercial insurance options.



Truncal
Garments



AIROS 6-Chamber & 8-Chamber Trunca Garment Sizing Chart

Model (CODE)	Inseam (CM)	Abdomen Length (CM)	Total Length (CM)	Ankle Circ. (CM)	Knee Circ. (CM)	Thigh Circ. (CM)	Abdomen Circ. (CM)	Waist/Hip Circ. (CM)	Belly Circ. (CM)
6-Chamber Small Pants (6-SP)	66	46	117	53	79	79	123	127.4	125.3
6-Chamber Medium Pants (6-MP)	81	46	127	64	81	81	125	127.6	123.8
6-Chamber Large Pants (6-LP)	91	46	142	71	94	94	145	148.3	145.2
8-Chamber Small Pants (8-SP)	75	41.9	116.9	54	76.4	79	122.8	127.4	125.3
8-Chamber Small Pants XL (8-SPXL)	75.4	41.5	116.9	53	79	88	155.7	155	155.7
8-Chamber Small Pants XXL (8-SPXXL)	75.4	41.5	116.9	53	79	93	170.7	170	170.7
8-Chamber Medium Pants (8-MP)	83.8	42.5	126.3	65.2	80	82	122	127.6	123.8
8-Chamber Medium Pants XL (8-MPXL)	83.8	42.4	126.2	64	82	90	155.9	155	155.9
8-Chamber Medium Pants XXL (8-MPXXL)	83.8	42.4	126.2	64	82	95	171	170	171
8-Chamber Large Pants (8-LP)	94.7	47.4	142.1	71	93.1	95	143	148.3	145.2
8-Chamber Large Pants XL (8-LPXL)	96.1	46	142.1	71	94	102	156.4	155	156.4
8-Chamber Large Pants XXL (8-LPXXL)	96.1	46	142.1	71	94	110	171.2	170	171.1